


Anti-Jab1/Csn5/Cops5 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa, Chip - Q92905

Item Number: CM0007866



Introduction

Rabbit monoclonal antibody against JAB1/Csn5/Cops5 (COP9 signalosome complex subunit 5), validated for WB, IHC-P, IF/ICC, ELISA, and CHIP in human, mouse, and rat samples. Ideal for ubiquitin-proteasome pathway studies.

[Learn More](#)

Parameter	Value
Product Name	JAB1/Csn5/Cops5 Rabbit mAb
Remarks/Alias	CSN5; JAB1; SGN5; MOV-34; JAB1/Csn5/Cops5
Species	Human
Gene ID	10987
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human JAB1/Csn5/Cops5 (Q92905)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, CHIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92905
Protein Weight	38kDa
Shipping	Ice bag