

Anti-Jab1/Csn5/Cops5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92905

Item Number: CM0016340



Introduction

Rabbit polyclonal antibody against human JAB1/Csn5/Cops5, a COP9 signalosome subunit involved in deneddylation and protein degradation. Validated for WB, IF/ICC, ELISA with reactivity in human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	JAB1/Csn5/Cops5 Rabbit pAb
Remarks/Alias	CSN5; JAB1; SGN5; MOV-34; JAB1/Csn5/Cops5
Species	Human
GeneID (Human)	10987
GeneID	10987
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human JAB1/Csn5/Cops5 (NP_006828.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92905
Protein Weight	38kDa
Shipping	Ice bag