


Anti-Jade1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Ie81

Item Number: CM0027625



Introduction

Rabbit polyclonal antibody targeting human JADE1, a scaffold subunit of HBO1 histone acetyltransferase complex. Validated in WB, IF/ICC, ELISA; cross-reacts with mouse. Ideal for chromatin and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	JADE1 Rabbit pAb
Remarks/Alias	PHF17, JADE1
Species	Human
Gene ID (Human)	79960
Gene ID	79960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-509 of human JADE1 (NP_079176.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6IE81
Protein Weight	96kDa
Shipping	Ice bag