


Anti-Jade1 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q6le81

Item Number: CM0017394



Introduction

JADE1 rabbit polyclonal antibody validated for IF/ICC and ELISA. Detects human and mouse JADE1, a scaffold protein of HBO1 complex involved in histone acetylation, apoptosis, and Wnt signaling. Ideal for chromatin research.

[Learn More](#)

Parameter	Value
Product Name	JADE1 Rabbit pAb
Remarks/Alias	PHF17; JADE1
Species	Human
Gene ID (Human)	79960
Gene ID	79960
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 210-509 of human JADE1 (NP_079176.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6IE81
Protein Weight	96kDa
Shipping	Ice bag