

Anti-Jag2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y219

Item Number: CM0018458



Introduction

JAG2 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets Protein jagged-2, a Notch ligand involved in development and signaling. Ideal for developmental biology and muscular dystrophy research.

[Learn More](#)

Parameter	Value
Product Name	JAG2 Rabbit pAb
Remarks/Alias	HJ2; SER2; LGMDR27; JAG2
Species	Human
Gene ID (Human)	3714
Gene ID	3714
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 24-200 of human JAG2 (NP_002217.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y219
Protein Weight	133kDa
Shipping	Ice bag