


Anti-Jak1 Polyclonal Antibody For Wb, If/Icc, Elisa - P23458

Item Number: CM0018794



Introduction

JAK1 Rabbit Polyclonal Antibody validated for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat JAK1. Suitable for cytokine signaling and JAK-STAT pathway studies.

[Learn More](#)

Parameter	Value
Product Name	JAK1 Rabbit pAb
Remarks/Alias	JTK3, AIIDE, JAK1A, JAK1B, JAK1
Species	Human
GeneID	3716
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 30-300 of human JAK1 (NP_002218.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23458
Protein Weight	133kDa
Shipping	Ice bag