

Anti-Jak1 Rabbit Monoclonal Antibody For Wb And Elisa - P23458

Item Number: CM0004638

Introduction

Rabbit monoclonal antibody against human JAK1, suitable for WB and ELISA. Cross-reacts with human, mouse, and rat. JAK1 is a tyrosine-protein kinase involved in cytokine signaling (P23458).

[Learn More](#)

Parameter	Value
Product Name	JAK1 Rabbit mAb
Remarks/Alias	JTK3; AIIDE; JAK1A; JAK1B
Species	Human
GeneID	3716
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 721-820 of human JAK1 (NP_002218.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23458
Protein Weight	133kDa
Shipping	Ice bag