

Anti-Jak1 Rabbit Monoclonal Antibody For Wb And Elisa - P23458

Item Number: CM0003746



Introduction

High-quality rabbit monoclonal antibody against human JAK1 tyrosine kinase (P23458). Validated for Western blot and ELISA; cross-reacts with mouse and rat. Ideal for cytokine signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	JAK1 Rabbit mAb
Remarks/Alias	JTK3; AIIDE; JAK1A; JAK1B; JAK1
Species	Human
GeneID (Human)	3716
GeneID	3716
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 600-700 of human JAK1 (P23458).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23458
Protein Weight	133kDa
Shipping	Ice bag