


Anti-Jak3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P52333

Item Number: CM0014086



Introduction

High-quality rabbit polyclonal antibody targeting human JAK3, a non-receptor tyrosine kinase involved in cytokine signaling. Validated in WB, IF/ICC, and ELISA. Suitable for immunology research.

[Learn More](#)

Parameter	Value
Product Name	JAK3 Rabbit pAb
Remarks/Alias	JAKL; L-JAK; JAK-3; L-JAK; JAK3_HUMAN; JAK3
Species	Human
Gene ID (Human)	3718
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-530 of JAK3 (NP_000206.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P52333
Protein Weight	125kDa
Shipping	Ice bag