

Anti-Jakmip2 Polyclonal Antibody For Wb, Elisa - Q96Aa8

Item Number: CM0015250

Introduction

Rabbit polyclonal antibody targeting JAKMIP2, validated for WB and ELISA. Reacts with human, mouse, and rat. Suitable for studying Golgi-associated coiled-coil protein in neuroendocrine tissues.

[Learn More](#)

Parameter	Details
Product Name	JAKMIP2 Rabbit pAb
Remarks/Alias	NECC1; JAMIP2; JAKMIP2
Species	Human
Gene ID (Human)	9832
Gene ID	9832
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 561-810 of human JAKMIP2 (NP_055605.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AA8
Protein Weight	95kDa
Shipping	Ice bag