

Anti-Jam-A/Cd321 Rabbit Polyclonal Antibody For Wb, Elisa - Q92993

Item Number: CM0018909



Introduction

Rabbit polyclonal antibody targeting JAM-A/CD321, suitable for WB and ELISA applications. Reacts with human, mouse, and rat samples. Useful for cell adhesion and tight junction studies. SWISS: Q92993.

[Learn More](#)

Parameter	Value
Product Name	JAM-A/CD321 Rabbit pAb
Remarks/Alias	JAM; Jcam; JAM-1; JAM-A; Jcam1; Ly106; ESTM33; 9130004G24; JAM-A/CD321
Species	Mouse
GeneID (human)	16456
GeneID	16456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-242 of mouse JAM-A/CD321. (NP_001193762.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92993
Protein Weight	32kDa
Shipping	Ice bag