

Anti-Jam-A/Cd321/F11R Monoclonal Antibody For Fc - Q9Y624

Item Number: CM0004178



Introduction

Rabbit monoclonal antibody targeting human JAM-A/CD321/F11R, suitable for flow cytometry. Recombinant immunogen (aa 28-242), cross-reacts with human. Ideal for tight junction, transmigration, and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human JAM-A/CD321/F11R mAb
Remarks/Alias	JAM; KAT; JAM1; JAMA; JCAM; CD321; PAM-1
Species	Human
Gene ID (Human)	50848
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 28-242 of human JAM-A/CD321/F11R (NP_058642.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9Y624
Protein Weight	27kDa/32kDa
Shipping	Ice bag