

Anti-Jam-A/Cd321/F11R Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Y624

Item Number: CM0007674



Introduction

Recombinant rabbit monoclonal antibody targeting human JAM-A (F11R, CD321). Validated for WB, IF/ICC, and ELISA. Specific for human. 33 kDa. Ideal for tight junction research, leukocyte transmigration studies, and viral entry analysis.

[Learn More](#)

Parameter	Value
Product Name	JAM-A/CD321/F11R Rabbit mAb
Remarks/Alias	JAM, KAT, JAM1, JAMA, JCAM, CD321, PAM-1, JAM-A/CD321/F11R
Species	Human
GeneID (Human)	50848
Immunogen	Recombinant protein of human JAM-A/CD321/F11R
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y624
Protein Weight	33kDa
Shipping	Ice bag