

Anti-Jam-A/Cd321/F11R Rabbit Polyclonal Antibody For Wb, IHC-P, Elisa - Q9Y624

Item Number: CM0014384



Introduction

Rabbit polyclonal antibody against human JAM-A (Junctional adhesion molecule A) for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Useful for tight junction and immunology research.

[Learn More](#)

Parameter	Value
Product Name	JAM-A/CD321/F11R Rabbit pAb
Remarks/Alias	JAM; KAT; JAM1; JAMA; JCAM; CD321; PAM-1; JAM-A/CD321/F11R
Species	Human
Gene ID	50848 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-238 of human JAM-A/CD321/F11R (NP_058642.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y624
Protein Weight	33kDa
Shipping	Ice bag