

Anti-Jam-A/Cd321/F11R Rabbit Polymab® Antibody For Wb, If/Icc, Elisa - Q9Y624 / O88792

Item Number: CM0005712



Introduction

Rabbit PolymAb® targeting JAM-A/CD321/F11R, validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat JAM-A. Ideal for research on tight junctions, platelet activation, and viral entry. UniProt: Q9Y624 / O88792.

[Learn More](#)

Parameter	Value
Product Name	JAM-A/CD321/F11R Rabbit PolymAb®
Remarks/Alias	JAM; Jcam; JAM-1; JAM-A; Jcam1; Ly106; ESTM33; 9130004G24; JAM-A/CD321
Species	Human/Mouse
Gene ID (Human)	50848
Gene ID	50848/16456
Immunogen	Recombinant protein of human/mouse JAM-A/CD321/F11R
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y624/O88792
Protein Weight	32kDa
Shipping	Ice bag