


## Anti-Jarid1B Polyclonal Antibody For Wb, Elisa - Q9Ugl1

Item Number: CM0015025



### Introduction

JARID1B rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat JARID1B. Ideal for chromatin biology and cancer research. - Q9UGL1

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | JARID1B Rabbit pAb                                                                                                     |
| Remarks/Alias    | CT31; PLU1; PUT1; MRT65; PLU-1; JARID1B; PPP1R98; RBP2-H1; RBBP2H1A                                                    |
| Species          | Human                                                                                                                  |
| Gene ID (Human)  | 10765                                                                                                                  |
| Gene ID          | 10765                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1400-1500 of human JARID1B (NP_006609.3) |
| Source           | Rabbit                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                  |
| Application      | WB, ELISA                                                                                                              |
| Cross Reactivity | Human, Mouse, Rat                                                                                                      |
| SWISS            | Q9UGL1                                                                                                                 |
| Protein Weight   | 176kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                                |