

Anti-Jarid2 Rabbit Monoclonal Antibody For Wb, Elisa - Q92833

Item Number: CM0008342



Introduction

Rabbit monoclonal antibody targeting JARID2 (Jumonji/ARID domain-containing protein 2), suitable for WB and ELISA. Cross-reacts with human, mouse, rat. Recognizes a PRC2 accessory subunit critical for embryonic development and chromatin regulation.

[Learn More](#)

Parameter	Value
Product Name	JARID2 Rabbit mAb
Remarks/Alias	JMJ; DIDDF; JARID2
Species	Human
GeneID (Human)	3720
GeneID	3720
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human JARID2 (Q92833).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92833
Protein Weight	139kDa
Shipping	Ice bag