

Anti-Jdp2 Polyclonal Antibody For Wb And Elisa - Q8Wyk2

Item Number: CM0014403

Introduction

Rabbit polyclonal antibody against human JDP2, validated for WB and ELISA, cross-reacting with rat. Suitable for studying AP-1 transcriptional repression and chromatin modification mechanisms.

[Learn More](#)

Parameter	Value
Product Name	JDP2 Rabbit pAb
Remarks/Alias	JUNDM2; JDP2
Species	Human
Gene ID	122953
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-174 of human JDP2 (NP_001128521.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8WYK2
Protein Weight	19kDa
Shipping	Ice bag