

Anti-Jmjd4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H9V9

Item Number: CM0029925

Introduction

Rabbit polyclonal antibody against human JMJD4, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, rat. Target involved in translation termination through ETF1 hydroxylation.

[Learn More](#)

Parameter	Value
Product Name	JMJD4 Rabbit pAb
Remarks/Alias	JMJD4
Species	Human
Gene ID (Human)	65094
Gene ID	65094
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-463 of human JMJD4 (NP_075383.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9V9
Protein Weight	52kDa
Shipping	Ice bag