

Anti-Jmjd6 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Nyc1

Item Number: CM0021524



Introduction

Rabbit polyclonal antibody against human JMJD6 (PSR/PTDSR), validated for WB, IHC-P, and ELISA. Reacts with human and mouse JMJD6, a bifunctional arginine demethylase and lysyl-hydroxylase involved in RNA splicing and transcriptional regulation.

[Learn More](#)

Parameter	Description
Product Name	JMJD6 Rabbit pAb
Remarks/Alias	PSR; PTDSR; PTDSR1; JMJD6
Species	Human
GeneID (Human)	23210
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human JMJD6 (NP_055982.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6NYC1
Protein Weight	46kDa
Shipping	Ice bag