

Anti-Jmjd7 Polyclonal Antibody For Wb, Elisa - P0C870

Item Number: CM0028549



Introduction

High-quality rabbit polyclonal antibody against human JMJD7 (P0C870), validated for WB and ELISA. Recognizes human, mouse, and rat JMJD7, a bifunctional peptidase and lysyl hydroxylase involved in histone cleavage and translation regulation.

[Learn More](#)

Parameter	Specification
Product Name	JMJD7 Rabbit pAb
Remarks/Alias	JMJD7
Species	Human
GeneID (Human)	100137047
GeneID	100137047
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-316 of human JMJD7 (NP_001108104.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0C870
Protein Weight	36kDa
Shipping	Ice bag