

Anti-Jmjd8 Polyclonal Antibody For Wb, Elisa - Q96S16

Item Number: CM0015961

Introduction

Rabbit polyclonal antibody targeting JMJD8 (JmjC domain-containing protein 8) for WB and ELISA applications. Reacts with human, mouse, rat. Involved in NF-kB signaling and cellular metabolism.

[Learn More](#)

Parameter	Value
Product Name	JMJD8 Rabbit pAb
Remarks/Alias	PP14397; C16orf20; JMJD8
Species	Human
GeneID (Human)	339123
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human JMJD8 (NP_001005920.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96S16
Protein Weight	30kDa
Shipping	Ice bag