

Anti-Jnk1 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P45983

Item Number: CM0006744



Introduction

Rabbit monoclonal antibody against human JNK1 (MAPK8, SAPK1c), suitable for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. For studying stress-activated kinase signaling, apoptosis, and cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	JNK1 Rabbit mAb
Remarks/Alias	JNK; JNK1; PRKM8; SAPK1; JNK-46; JNK1A2; SAPK1c; JNK21B1/2
Species	Human
Gene ID (Human)	5599
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human JNK1 (NP_620637.1)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45983
Protein Weight	48kDa
Shipping	Ice bag