

## Anti-Jnk1 Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ip, Elisa - P45983

Item Number: CM0007124



### Introduction

Anti-JNK1 rabbit monoclonal antibody [KD Validated] for WB, IP, ELISA. Detects human, mouse, and rat JNK1 (SAPK1c, MAPK8) with a molecular weight of 48 kDa. Validated by knock-down, suitable for signaling research.

[Learn More](#)

| Parameter        | Value                                                                                           |
|------------------|-------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] JNK1 Rabbit mAb                                                                  |
| Remarks/Alias    | JNK; [KD Validated] JNK1                                                                        |
| Species          | Human                                                                                           |
| GeneID           | 5599                                                                                            |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 250-350 of human JNK1 (P45983) |
| Source           | Rabbit                                                                                          |
| Category         | Monoclonal Antibodies                                                                           |
| Application      | WB, IP, ELISA                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                               |
| SWISS            | P45983                                                                                          |
| Protein Weight   | 48kDa                                                                                           |
| Shipping         | Ice bag                                                                                         |