

Anti-Jnk1 Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - P45983

Item Number: CM0016304



Introduction

Rabbit polyclonal JNK1 antibody validated for WB, IF/ICC, IHC-P, and ELISA. Detects human, mouse, rat MAPK8/SAPK1c kinase involved in stress signaling, apoptosis, and autophagy. Ideal IVD assay development tool.

[Learn More](#)

Parameter	Value
Product Name	JNK1 Rabbit pAb
Remarks/Alias	JNK; JNK1; PRKM8; SAPK1; JNK-46; JNK1A2; SAPK1c; JNK21B1/2
Species	Human
GeneID (Human)	5599
GeneID	5599
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 245-345 of human JNK1 (NP_620635.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45983
Protein Weight	48 kDa
Shipping	Ice bag