

Anti-Jnk1/2/3 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P45983 / P45984 / P53779

Item Number: CM0004038



Introduction

Rabbit monoclonal antibody targeting JNK1/2/3 (P45983/P45984/P53779) validated for Western blot, immunofluorescence, and ELISA in human, mouse, and rat samples. Ideal for stress signaling and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	JNK1/2/3 Rabbit mAb
Remarks/Alias	JNK; JNK-46; JNK1; JNK1A2; JNK21B1/2; PRKM8; SAPK1; SAPK1c; JNK1/2/3
Species	Human
GeneID	5599/5601/5602
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-383 of human JNK1 (P45983).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45983/P45984/P53779
Protein Weight	46,54kDa
Shipping	Ice bag