

Anti-Jnk1/2/3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P45983 / P45984 / P53779

Item Number: CM0006930



Introduction

Recombinant rabbit monoclonal antibody against JNK1/2/3 (MAPK8/9/10). Suitable for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets stress-activated kinases involved in apoptosis, proliferation, and circadian regulation. Ideal for signaling research.

[Learn More](#)

Parameter	Value
Product Name	JNK1/2/3 Rabbit mAb
Remarks/Alias	JNK; JNK1; PRKM8; SAPK1; JNK-46; JNK1A2; SAPK1c; JNK21B1/2; JNK3; JNK3A; PRKM10; SAPK1b; p493F12; p54bSAPK; JNK1/JNK3
Species	Human
GeneID (Human)	5599/5601/5602
GeneID	5599/5601/5602
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human JNK1 (NP_620637.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45983/P45984/P53779
Protein Weight	35kDa/44kDa/48kDa
Shipping	Ice bag