

Anti-Jnk1/2/3 Rabbit Polyclonal Antibody For Wb - P45983/P45984/P53779

Item Number: CM0020148



Introduction

Rabbit polyclonal antibody against JNK1/2/3 (MAPK8/9/10). Validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with Human, Mouse, Rat. Suitable for studying stress-activated kinase signaling, apoptosis, and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	JNK1/2/3 Rabbit pAb
Species	Human
Gene ID (Human)	5599/5601/5602
Gene ID	5599/5601/5602
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-370 of human JNK1/2/3 (NP_620637.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45983/P45984/P53779
Protein Weight	35kDa/44kDa/48kDa/27kDa/52kDa
Shipping	Ice bag