

Anti-Jnk1/Jnk3 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P45983 / P53779

Item Number: CM0004063



Introduction

High-quality rabbit monoclonal antibody against JNK1/JNK3 (P45983/P53779), validated for WB, IHC-P, IF/ICC, ELISA. Suitable for stress signaling, apoptosis, and inflammation research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	JNK1/JNK3 Rabbit mAb
Remarks/Alias	JNK; JNK-46; JNK1; JNK1A2; JNK21B1/2; PRKM8; SAPK1; SAPK1c; JNK1/JNK3
Species	Human
Gene ID (Human)	5599/5602
Gene ID	5599/5602
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-383 of human JNK1 (P45983). This antibody does not cross react with JNK2 protein.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45983/P53779
Protein Weight	35kDa/44kDa/48kDa
Shipping	Ice bag