


Anti-Jnk2 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - P45984

Item Number: CM0014858



Introduction

Rabbit polyclonal antibody against human JNK2 (MAPK9), validated for WB and ELISA. Detects JNK2 in human and mouse samples. Ideal for studying SAPK/JNK signaling, apoptosis, and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] JNK2 Rabbit pAb
Remarks/Alias	JNK2; SAPK; p54a; JNK2A; JNK2B; PRKM9; JNK-55; SAPK1a; JNK2BETA; p54aSAPK; JNK2ALPHA; K2
Species	Human
GeneID (Human)	5601
GeneID	5601
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P45984
Protein Weight	48kDa
Shipping	Ice bag