


Anti-Jph1 Rabbit Polyclonal Antibody For Wb - Q9Hdc5

Item Number: CM0020297



Introduction

High-quality rabbit polyclonal antibody against human JPH1 (Junctophilin-1), validated for Western blot and ELISA. Cross-reacts with mouse; ideal for skeletal muscle junctional complex studies.

[Learn More](#)

Parameter	Value
Product Name	JPH1 Rabbit pAb
Remarks/Alias	JP1; JP-1; CMT2K; JPH1
Species	Human
Gene ID (Human)	56704
Gene ID	56704
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 410-640 of human JPH1 (NP_065698.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9HDC5
Protein Weight	72kDa
Shipping	Ice bag