

Anti-Jph4 Polyclonal Antibody For Wb, Elisa - Q96Jj6

Item Number: CM0016070

Introduction

Rabbit polyclonal antibody against human Junctophilin-4 (JPH4). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for studying junctional membrane complex formation and neuronal calcium signaling.

[Learn More](#)

Parameter	Value
Product Name	JPH4 Rabbit pAb
Remarks/Alias	JP4; JPHL1; JPH4
Species	Human
GeneID (human)	84502
GeneID	84502
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 405-605 of human JPH4 (NP_115828.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Jj6
Protein Weight	66kDa
Shipping	Ice bag