

Anti-Jrkl Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y4A0

Item Number: CM0013126

Introduction

Rabbit polyclonal antibody against human JRKL (Jerky protein homolog-like), validated for WB and ELISA. Recognizes human and mouse JRKL. Ideal for nuclear protein studies and DNA-binding research.

[Learn More](#)

Parameter	Value
Product Name	JRKL Rabbit pAb
Remarks/Alias	HHMJG; JRKL
Species	Human
Gene ID (Human)	8690
Gene ID	8690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 225-524 of human JRKL (NP_003763.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y4A0
Protein Weight	60kDa
Shipping	Ice bag