

Anti-Junb Monoclonal Antibody For Wb, Ihc-P, Elisa - P17275

Item Number: CM0010949

Introduction

High-quality anti-JunB rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Targets human JunB transcription factor (AP-1 subunit); 36 kDa. Ideal for gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	JunB Rabbit mAb
Remarks/Alias	AP-1; JunB
Species	Human
Gene ID (Human)	3726
Gene ID	3726
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human JunB (P17275)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P17275
Protein Weight	36kDa
Shipping	Ice bag