

Anti-Junb Polyclonal Antibody For Wb, If/Icc, Elisa - P17275

Item Number: CM0027475

Introduction

Rabbit polyclonal antibody targeting human JunB (AP-1 transcription factor subunit) validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying gene regulation and cell signaling.

[Learn More](#)

Parameter	Value
Product Name	JunB Rabbit pAb
Remarks/Alias	AP-1; JunB
Species	Human
GeneID (Human)	3726
GeneID	3726
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human JunB (NP_002220.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17275
Protein Weight	36kDa
Shipping	Ice bag