

## Anti-Junb Rabbit Polyclonal Antibody For Wb, Elisa - P17275

Item Number: CM0017891

### Introduction

Rabbit polyclonal antibody targeting human JunB transcription factor, a key AP-1 component. Raised with recombinant immunogen; suitable for Western blot and ELISA applications. Ideal for studying gene regulation and transcription factor activity.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | JunB Rabbit pAb                                                                                                 |
| Remarks/Alias    | AP-1; JunB                                                                                                      |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 3726                                                                                                            |
| GeneID           | 3726                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human JunB (NP_002220.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, ELISA                                                                                                       |
| Cross Reactivity | Human                                                                                                           |
| SWISS            | P17275                                                                                                          |
| Protein Weight   | 36kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |