


Anti-JunD Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P17535

Item Number: CM0012225



Introduction

Recombinant anti-JunD rabbit monoclonal antibody validated for WB, IHC-P, IP, ELISA. Detects human, mouse, rat JunD (35kDa). For AP-1 transcription factor research.

[Learn More](#)

Parameter	Value
Product Name	JunD Rabbit mAb
Remarks/Alias	AP-1; JunD
Species	Human
GeneID (Human)	3727
GeneID	3727
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-141 of human JunD (P17535).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17535
Protein Weight	35kDa
Shipping	Ice bag