

Anti-K-Cadherin (Cdh6) Rabbit Monoclonal Antibody For Wb, Elisa - P55285

Item Number: CM0012567



Introduction

High-quality rabbit monoclonal antibody targeting human K-Cadherin (CDH6). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for cell adhesion, kidney, and brain research studies.

[Learn More](#)

Parameter	Value
Product Name	K-Cadherin (CDH6) Rabbit mAb
Remarks/Alias	CAD6; KCAD; K-Cadherin (CDH6)
Species	Human
GeneID	1004
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 691-790 of human K-Cadherin (CDH6) (P55285)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55285
Protein Weight	88kDa
Shipping	Ice bag