

## Anti-K-Cadherin (Cdh6) Rabbit Polyclonal Antibody For Wb, Elisa - P55285

Item Number: CM0029363



### Introduction

Rabbit polyclonal antibody targeting human K-Cadherin (CDH6), ideal for WB and ELISA. Cross-reacts with mouse and rat. Detects an 88kDa band. Suitable for cell adhesion and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                                      |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | K-Cadherin (CDH6) Rabbit pAb                                                                                                               |
| Remarks/Alias    | CAD6; KCAD; K-Cadherin (CDH6)                                                                                                              |
| Species          | Human                                                                                                                                      |
| GeneID (human)   | 1004                                                                                                                                       |
| GeneID           | 1004                                                                                                                                       |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human K-Cadherin (K-Cadherin (CDH6)) (NP_004923.1). |
| Source           | Rabbit                                                                                                                                     |
| Category         | Polyclonal Antibodies                                                                                                                      |
| Application      | WB, ELISA                                                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                          |
| SWISS            | P55285                                                                                                                                     |
| Protein Weight   | 88kDa                                                                                                                                      |
| Shipping         | Ice bag                                                                                                                                    |