

Anti-Kaiso/Zbtb33 Polyclonal Antibody - Q86T24

Item Number: CM0013696



Introduction

Kaiso/ZBTB33 Rabbit Polyclonal Antibody, validated for WB, ELISA, and ChIP applications. Detects endogenous Kaiso protein in human, mouse, and rat. Ideal for chromatin immunoprecipitation and gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Kaiso/ZBTB33 Rabbit pAb
Remarks/Alias	ZNF348; ZNF-kaiso; Kaiso/ZBTB33
Species	Human
Gene ID (Human)	10009
Gene ID	10009
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 523-672 of human Kaiso/Kaiso/ZBTB33 (NP_006768.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86T24
Protein Weight	74kDa
Shipping	Ice bag