

Anti-Kaiso/Zbtb33 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q86T24

Item Number: CM0012309



Introduction

Rabbit monoclonal antibody against human Kaiso/ZBTB33 (ZNF348), validated for WB, IF/ICC, IP, and ELISA. Reacts with human and rat samples. Suitable for studying transcriptional regulation, Wnt signaling, and epigenetic mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Kaiso/ZBTB33 Rabbit mAb
Remarks/Alias	ZNF348; ZNF-kaiso; Kaiso/ZBTB33
Species	Human
Gene ID (Human)	10009
Gene ID	10009
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 573-672 of human Kaiso/ZBTB33 (Q86T24).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	Q86T24
Protein Weight	74KDa
Shipping	Ice bag