

## Anti-Kal1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P23352

Item Number: CM0019956



### Introduction

High-quality rabbit polyclonal antibody against human Anosmin-1 (KAL1). Validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for studying olfactory bulb development and Kallmann syndrome.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | KAL1 Rabbit pAb                                                                                                           |
| Remarks/Alias    | HH1; HHA; KAL; KMS; KAL1; ADMLX; WFDC19; KALIG-1                                                                          |
| Species          | Human                                                                                                                     |
| Gene ID (Human)  | 3730                                                                                                                      |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human KAL1 (NP_000207.2). |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, IF-P, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | P23352                                                                                                                    |
| Protein Weight   | 76kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |