


Anti-Kallikrein 10 (Nes1/Prssl1) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O43240

Item Number: CM0009525



Introduction

High-specificity rabbit monoclonal antibody targeting human Kallikrein 10 (NES1/PRSSL1; O43240). Suitable for Western blot, IF/ICC, and ELISA. Ideal for breast, ovarian, and prostate cancer biomarker research.

[Learn More](#)

Parameter	Value
Product Name	Kallikrein 10/K10 Rabbit mAb
Remarks/Alias	NES1; PRSSL1
Species	Human
Gene ID (Human)	5655
Gene ID	5655
Immunogen	Recombinant protein corresponding to amino acids 43-276 of human KLK10 (NP_665895.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O43240
Protein Weight	30 kDa
Shipping	Ice bag