


Anti-Kallikrein-5 (Klk5) Polyclonal Antibody For Wb And Elisa - Q9Y337

Item Number: CM0027709



Introduction

High-quality rabbit anti-KLK5 polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for desquamation studies and serine protease research.

[Learn More](#)

Parameter	Value
Product Name	KLK5 Rabbit pAb
Remarks/Alias	SCTE; KLKL2; KLK-L2; KLK5
Species	Human
GeneID (Human)	25818
GeneID	25818
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 41-177 of human KLK5 (NP_036559.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y337
Protein Weight	32kDa
Shipping	Ice bag