

Anti-Kallikrein 8 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O60259

Item Number: CM0005914



Introduction

Rabbit monoclonal antibody against human Kallikrein 8 (Neurotrophin) validated for Western blot, IHC-P, and ELISA. Suitable for neuroscience, skin biology, and cancer research. Reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	Kallikrein 8 Rabbit mAb
Remarks/Alias	NP; HNP; NRPN; PRSS19; TADG14; Kallikrein 8
Species	Human
GeneID (human)	11202
GeneID	11202
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 161-260 of human Kallikrein 8 (O60259).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O60259
Protein Weight	28kDa
Shipping	Ice bag