

Anti-Kank2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q63Zy3

Item Number: CM0021372



Introduction

Rabbit polyclonal antibody against human KANK2, suitable for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Targets a protein involved in transcription regulation, apoptosis, and actin cytoskeleton organization.

[Learn More](#)

Parameter	Value
Product Name	KANK2 Rabbit pAb
Remarks/Alias	SIP; MXRA3; PPKWH; NPHS16; ANKRD25; KANK2
Species	Human
Gene ID (Human)	25959
Gene ID	25959
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human KANK2 (NP_001129663.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q63ZY3
Protein Weight	91kDa
Shipping	Ice bag