

Anti-Kansl3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9P2N6

Item Number: CM0030889

Introduction

Rabbit polyclonal antibody targeting human KANSL3, validated for Western blot and ELISA. KANSL3 is a key component of the NSL histone acetyltransferase complex, involved in chromatin regulation and microtubule stabilization. Ideal for epigenetics and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	KANSL3 Rabbit pAb
Remarks/Alias	NSL3; Rcd1; KIAA1310; KANSL3
Species	Human
GeneID (human)	55683
GeneID	55683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human KANSL3 (NP_001108488.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9P2N6
Protein Weight	96kDa
Shipping	Ice bag