


Anti-Kap1/Trim28 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13263

Item Number: CM0002591



Introduction

Anti-KAP1/TRIM28 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Recognizes the transcription corepressor and E3 ligase TRIM28; molecular weight ~89 kDa.

[Learn More](#)

Parameter	Value
Product Name	KAP1/TRIM28 Rabbit mAb
Remarks/Alias	KAP1; TF1B; RNF96; TIF1B; PPP1R157; TIF1beta; KAP1/TRIM28
Species	Human
GeneID (human)	10155
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human KAP1/TRIM28 (Q13263).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13263
Protein Weight	89 kDa
Shipping	Ice bag