


Anti-Kap1/Trim28 Polyclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q13263

Item Number: CM0020732



Introduction

High-quality anti-KAP1/TRIM28 rabbit polyclonal antibody validated for WB, IP, IF/ICC, IHC-P, and ELISA. Detects KAP1/TRIM28 in human, mouse, and rat. Ideal for chromatin biology, transcription regulation, and immune signaling research.

[Learn More](#)

Parameter	Specification
Product Name	KAP1/TRIM28 Rabbit pAb
Remarks/Alias	KAP1; TF1B; RNF96; TIF1B; PPP1R157; TIF1beta; KAP1/TRIM28
Species	Human
Gene ID	10155
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 576-835 of human KAP1/TRIM28 (NP_005753.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13263
Protein Weight	89kDa
Shipping	Ice bag