

Anti-Karyopherin B3 (Ipo5) Monoclonal Antibody For Wb, Elisa - O00410

Item Number: CM0007438



Introduction

Recombinant rabbit monoclonal antibody specific for human Karyopherin β 3 (IPO5/Importin-5). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for nuclear transport and importin research.

[Learn More](#)

Parameter	Value
Product Name	Karyopherin β 3 (IPO5) Rabbit mAb
Remarks/Alias	IMB3; Pse1; imp5; KPNB3; RANBP5; Karyopherin β 3 (IPO5)
Species	Human
GeneID	3843
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human Karyopherin β 3 (IPO5). (NP_002262.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00410
Protein Weight	124kDa
Shipping	Ice bag